




SENDING DATA TO YOUR AIR LIQUIDE CONSULTANT


Your device is equipped internally with bluetooth and 4G cellular capabilities. This enables your device to connect to cellular networks in your area and remotely **sends your therapy information/data seamlessly to your AirLiquide provider**. Although your device will automatically send data daily, you **Air Liquide consultant may ask you to send the data manually or bring your device to the clinic**.


Data transfer via iDream


- 1.** 

From the main screen, press the **menu key**.
- 2.** 

Press the **- key** to access the device information.
- 3.** 

Use the **menu key** to move left until the data "send" screen.
- 4.** 

Use the **+ key** to change 'no' to 'yes'.
- 5.** 

Finalise the process by using the **ramp key** to initiate data transfer.
- 6.** 

Allow time for the device to connect and transfer data over to the cloud

Data transfer via PrismaApp on your smartphone

Once your device and smartphone have been paired via bluetooth, open the application and click on the last tab on the bottom of the screen named **'other'**.

You can activate the manual data transfer by pressing **'Retrieve new therapy data from the device'**

